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

Fig. 1: Subgroups of patients with chronic heart failure.

<sup>1</sup> Facts and Figures: Heart Failure. European Society of Cardiology website. [http://www.escardio.org/static\\_file/Escardio/Press-media/Facts-Figures-HF09.pdf](http://www.escardio.org/static_file/Escardio/Press-media/Facts-Figures-HF09.pdf). Accessed 22 June 2015.

<sup>2</sup> Fauci AS, Braunwald E, Kasper DL, et al, eds. Harrison's Principles of Internal Medicine. 17th ed. New York: McGraw-Hill; 2008.

<sup>3</sup> Levy D, Kenchaiah S, Larson MG, et al. Long-term trends in the incidence of and survival with heart failure. *N Engl J Med*. 2002;347(18):1397-1402.

<sup>4</sup> Stewart S, MacIntyre K, Hole DJ, Capewell S, McMurray JJ. More 'malignant' than cancer? Five-year survival following a first admission for heart failure. *Eur J Heart Fail*. 2001; 3:315-322.

<sup>5</sup> Fauci AS, Braunwald E, Kasper DL, et al, eds. Harrison's Principles of Internal Medicine. 17th ed. New York: McGraw-Hill; 2008.

<sup>6</sup> Cowie MR, Anker SD, Cleland JGF. Improving Care for Patients With Acute Heart Failure: Before, During and After Hospitalization. Oxford PharmaGenesis; 2014. <http://www.oxfordhealthpolicyforum.org/AHFreport>. Accessed June 22, 2015.

<sup>7</sup> J.J.V. McMurray, et al. (2012) ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure 2012. *European Heart Journal* 2012; 33, 1787-1847.

<sup>8</sup> P. Ponikowski, et al. 2016 Guidelines for the diagnosis and treatment of acute and chronic heart failure. *European Heart J* 2016; 37, 2129-2200

<sup>9</sup> L.H. Lund, et al. (2012) Prevalence, correlates, and prognostic significance of QRS prolongation in heart failure with reduced and preserved ejection fraction. *European Heart J* 2013; 34: 529-539.

<sup>10</sup> K.-H. Kuck, et al. (2014) New devices in heart failure: an European Heart Rhythm Association report: developed by the European Heart Rhythm Association; endorsed by the Heart Failure Association. *Europace* 2014; 16: 109-128.

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

Fig. 2: CCM improves cardiac function and leads to reverse remodeling in patients.<sup>12</sup>

<sup>11</sup> W.T. Abraham, S.A. Smith (2013) Devices in the management of advanced, chronic heart failure. *Nat Rev Cardiol.* 2013; 10: 98–110.  
<sup>12</sup> M. Borggrefe, D. Burkhoff (2012) Clinical effects of cardiac contractility modulation (CCM) as a treatment for chronic heart failure. *Eur J Heart Fail* 2012; 14: 703–712.  
<sup>13</sup> C.M. Yu, et al. (2009) Impact of cardiac contractility modulation on left ventricular global and regional function and remodeling. *JACC Cardiovasc Imaging* 2009; 2: 1341–1349.  
<sup>14</sup> C. Butter, et al. (2008) Cardiac contractility modulation electrical signals improve myocardial gene expression in patients with heart failure. *J Am Coll Cardiol.* 2008;51(18):1784-9.  
<sup>15</sup> AR. Lyon, et al. (2013) Cardiac contractility modulation therapy in advanced systolic heart failure. *Nat Rev Cardiol.* 2013;10(10):584-98.  
<sup>16</sup> F. Giallauria, et al. (2014) Effects of cardiac contractility modulation by non-excitatory electrical stimulation on exercise capacity and quality of life: an individual patient's data meta-analysis of randomized controlled trials. *Int J Cardiol* 175 (2): 352–357.  
<sup>17</sup> J. Kuschyk, et al. (2015) Efficacy and survival in patients with cardiac contractility modulation: Long-term single center experience in 81 patients. *Int J Cardiol* 2015 (183C): 76–81.  
<sup>18</sup> A. Kloppe, et al. (2016) Long-term survival with Cardiac Contractility Modulation in patients with NYHA II or III symptoms and normal QRS duration. *Int J Cardiol.*2016; 21968.

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
Fig. 3: Potential therapy flow chart for this group of patients.

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<sup>20</sup> RC Gupta, et al. (2009) Ca(2+)-binding proteins in dogs with heart failure: effects of cardiac contractility modulation electrical signals. Clin Transl Sci. 2009; 2(3):211-5.

<sup>21</sup> C Butter, et al. (2008) Cardiac contractility modulation electrical signals improve myocardial gene expression in patients with heart failure. J Am Coll Cardiol. 2008; 51(18):1784-9.

<sup>22</sup> C Butter, et al. (2007) Enhanced inotropic state of the failing left ventricle by cardiac contractility modulation electrical signals is not associated with increased myocardial oxygen consumption. J Card Fail. 2007; 13(2):137-42.

<sup>23</sup> G Goliash, et al. (2012) The effect of device-based cardiac contractility modulation therapy on myocardial efficiency and oxidative metabolism in patients with heart failure. Eur J Nucl Med Mol Imaging. 2012; 39(3):408-15.

<sup>24</sup> RC Gupta, et al. (2009) Ca(2+)-binding proteins in dogs with heart failure: effects of cardiac contractility modulation electrical signals. Clin Transl Sci. 2009; 2(3):211-5.

<sup>25</sup> C Butter, et al. (2007) Enhanced inotropic state of the failing left ventricle by cardiac contractility modulation electrical signals is not associated with increased myocardial oxygen consumption. J Card Fail. 2007; 13(2):137-42.

<sup>26</sup> AR. Lyon, et al. (2013) Cardiac contractility modulation therapy in advanced systolic heart failure. Nat Rev Cardiol. 2013;10(10):584-98.

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


Fig. 4 and 5: Data from the FIX-HF-5 Trial. Change in maximum oxygen uptake ( $\Delta$  Peak  $VO_2$ ) and quality of life (MLWHFQ) respectively.

<sup>27</sup> A. Kadish, et al. (2011) A randomized controlled trial evaluating the safety and efficacy of cardiac contractility modulation in advanced heart failure. Am Heart J. 2011; 161(2):329-337.e1-2.

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
Fig. 6: CCM leads to greater improvement in function and quality of life.<sup>30, 31, 32</sup> (for CCM); <sup>33, 34, 35</sup> (for CRT)

<sup>28</sup> M. Borggreffe, D. Burkhoff (2012) Clinical effects of cardiac contractility modulation (CCM) as a treatment for chronic heart failure. Eur J Heart Fail 2012; 14: 703–712.

<sup>29</sup> W.T. Abraham, et al. (2011) Subgroup analysis of a randomized controlled trial evaluating the safety and efficacy of cardiac contractility modulation in advanced heart failure. J Card Fail. 2011; 17(9):710-7.

<sup>30</sup> W.T. Abraham, et al. (2011) Subgroup analysis of a randomized controlled trial evaluating the safety and efficacy of cardiac contractility modulation in advanced heart failure. J Card Fail. 2011; 17(9):710-7.

<sup>31</sup> D. Burkhoff, et al. (2010) Cardiac contractility modulation by electrical signals applied during the absolute refractory period as a treatment for chronic heart failure. Heart Failure Device Management. Edited by AM Feldman. 2010;44-58.

<sup>32</sup> M. Borggreffe, D. Burkhoff (2012) Clinical effects of cardiac contractility modulation (CCM) as a treatment for chronic heart failure. Eur J Heart Fail 2012; 14: 703–712.

<sup>33</sup> W.T. Abraham, et al. (2002) Cardiac Resynchronization in Chronic Heart Failure. N Engl J Med 2002; 346:1845-1853.

<sup>34</sup> J.B. Young, et al. (2003) Combined cardiac resynchronization and implantable cardioversion defibrillation in advanced chronic heart failure: the MIRACLE ICD Trial. JAMA 2003; 289(20): 2685-94.

<sup>35</sup> S. Chen, et al. (2012) Effect of cardiac resynchronization therapy and implantable cardioverter defibrillator on quality of life in patients with heart failure: a meta-analysis. Europace. 2012 Nov;14(11):1602-7.

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
Fig. 7: Giallauria et al. (IJC 2014) Effects of cardiac contractility modulation by non-excitatory electrical stimulation on exercise capacity and quality of life: an individual patient's data meta-analysis of randomized controlled trials.

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
Fig. 8: Observed Kaplan-Meier survival curve (blue) compared to point estimates of survival at 1 and 3 years provided per patient by the MAGGIC score.

<sup>36</sup> J. Kuschyk, et al. (2015) Efficacy and survival in patients with cardiac contractility modulation: Long-term single center experience in 81 patients. Int J Cardiol 2015 (183C): 76–81.

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Table 1. Efficacy parameters at baseline and at long-term follow-up (mean and SD).<sup>38</sup>

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NT-proBNP: N-terminales Propeptid BNP (Brain Natriuertic Peptide)

## Beachtenswerte Publikationen

- **Abi-Samra F, Gutterman D. Cardiac contractility modulation: a novel approach for the treatment of heart failure.** *Heart Failure Review*, 2016
- **Ponikowski P, A. Voors, S. Anker, et al. 2016 Guidelines for the diagnosis and treatment of acute and chronic heart failure.** *European Heart Journal*, 2016
- **Kloppe A, Lawo T, Mijic D, et al. Long-term survival with Cardiac Contractility Modulation in patients with NYHA II or III symptoms and normal QRS duration.** *International Journal of Cardiology*, 2016
- **Liu Ming, Fang Fang, Luo Xiu Xia, et al. Improvement of long-term survival by cardiac contractility modulation in heart failure patients: A case-control study.** *International Journal of Cardiology*, 2016
- **Kuschyk J, Roeger S, Schneider R, et al. Efficacy and survival in patients with cardiac contractility modulation: Long-term single center experience in 81 patients.** *International Journal of Cardiology*, 2015
- **Giallauria F, Vigorito C, Piepoli MF, Stewart Coats AJ. Effects of cardiac contractility modulation by non-excitatory electrical stimulation on exercise capacity and quality of life: an individual patient's data meta-analysis of randomized controlled trials.** *International Journal of Cardiology*, 2014
- **Kuck KH, Bordachar P, Borggreffe M, et al. New devices in heart failure: an European Heart Rhythm Association report: developed by the European Heart Rhythm Association; endorsed by the Heart Failure Association.** *Europace*, 2014
- **Lyon AR, Samara MA, Feldman DS. Cardiac contractility modulation therapy in advanced systolic heart failure.** *Nature Reviews Cardiology*, 2013
- **Borggreffe M, Burkhoff D. Clinical effects of cardiac contractility modulation (CCM) as a treatment for chronic heart failure.** *European Journal of Heart Failure*, 2012
- **Abraham WT, Nademanee K, Volosin K, et al. Subgroup analysis of a randomized controlled trial evaluating the safety and efficacy of cardiac contractility modulation in advanced heart failure.** *Journal of Cardiac Failure*, 2011
- **Yu CM, Chan JY, Zhang Q, et al. Impact of cardiac contractility modulation on left ventricular global and regional function and remodeling.** *JACC Cardiovascular Imaging*, 2009
- **Butter C, Rastogi S, Minden HH, Meyhöfer J, Burkhoff D, Sabbah HN. Cardiac contractility modulation electrical signals improve myocardial gene expression in patients with heart failure.** *Journal of the American College of Cardiology*, 2008
- **Borggreffe MM, Lawo T, Butter C, et al. Randomized, double blind study of non-excitatory, cardiac contractility modulation electrical impulses for symptomatic heart failure.** *European Heart Journal*, 2008

<sup>37</sup> A. Kloppe, et al. (2016) Long-term survival with Cardiac Contractility Modulation in patients with NYHA II or III symptoms and normal QRS duration. *Int J Cardiol*.2016; 21968.

<sup>38</sup> J. Kuschyk, et al. (2015) Efficacy and survival in patients with cardiac contractility modulation: Long-term single center experience in 81 patients. *Int J Cardiol* 2015 (183C): 76-81.